

Produkt information

Comfort Focus



Improved concentration and learning for pupils with:

- Concentration problems
- Reading and writing difficulties
- APD (Auditory Processing Disorders)




Comfort Focus is a wireless digital system for daily use in a school environment for pupils with concentration problems and other learning difficulties. The system delivers a superior sound quality. This is possible thanks to Comfort Audio's unique technology which transfers the sound digitally between the teacher's microphone and the pupil's receiver.

Wearing a receiver the pupil finds it easier to concentrate and can focus on the teacher for a longer period.

In 2007, in a school in Norrköping, Sweden, a study was started how assistive listening devices might help children with cognitive disorders in school. In the chart below, pupils indicate how well they could concentrate with or without the wireless system.

In its report, the school expresses the following: "The teachers describe the pupils as happier and that they find it easier to concentrate. Another advantage is that they do not draw attention to the pupil by going over things again with them as much as before."





The chart shows that pupils, according to their own evaluation, find it easier to concentrate using the assistive listening device.



Comfort Audio
HEAR THE FUTURE

Product information (continued)

Comfort Focus



Teacher

The teacher uses a small, flexible microphone which is attached to his clothes with a clip. Alternatively, a clip-on microphone or a boom-microphone can be used. A teacher's unit can be used with any number of pupil receivers.



Pupil

The pupil/s wear their receivers around the neck on a neck loop and listen with head- or earphones. Alternatively, a receiver can be worn directly at the ear. The volume is easily adjusted by using buttons on the long side of the device. The maximum volume is restricted to avoid it being too loud.

The standard Comfort Focus package contains

- Teacher microphone and pupil receiver
- Two charging stands
- Two neck loops, one long, one short
- Earphones
- Headphones
- Bag
- Manual on USB memory drive
- Clip-on microphone

Available accessories/spare parts

- Extra pupil receivers
- Replacement cushions for headphones
- Boom-microphone
- Loudspeakers for the classroom