

# Product Overview

COMFORT DIGISYSTEM

## Microphone



Conference Microphone DC10



Conference Microphone DC20



Microphone DM05



Microphone DM10



Microphone DM20



Microphone DM30



Microphone DM90

## Receiver



Receiver DE10



Receiver DH10



Receiver DT20

### Features

Features	Conference Microphone DC10	Conference Microphone DC20	Microphone DM05	Microphone DM10	Microphone DM20	Microphone DM30	Microphone DM90	Receiver DE10	Receiver DH10	Receiver DT20
Type of microphone	Long range	Long range	Close range	Close range	Close range	Close range	Close range			
Integrated receiver		◆				◆	◆			
Eavesdrop-secure	◆	◆	◆	◆	◆	◆	◆	◆	◆	◆
Zoom-microphone	◆	◆					◆			
Microphone modes	4	4	1	1	1	1	2			
Transient control	◆	◆								
Noise reduction	◆	◆								
Display and language choice	◆	◆		◆	◆	◆	◆		◆	
Function indicator	◆	◆	◆	◆	◆	◆	◆	◆	◆	
Landline telephone		◆								
Mobile phone	◆	◆					◆			
Socket for headphones	◆	◆				◆	◆		◆	
Line-in	◆	◆		◆	◆	◆	◆			
Induction loop						◆	◆		◆	
Remote control		◆								
Volume control	◆	◆		◆	◆	◆	◆	◆	◆	
PIN-Code	◆	◆		◆	◆	◆	◆		◆	
Bluetooth	◆	◆					◆			
Automatic Gain Control (AGC)	◆	◆	◆	◆	◆	◆	◆			
Perceptual Speech Enhancement (PSE)	◆	◆								
SecureStream Technology (SST)	◆	◆	◆	◆	◆	◆	◆	◆	◆	◆
Free Frequency Check (FFC)	◆	◆		◆	◆	◆	◆		◆	
Free Channel Request (FCR)				◆	◆	◆				
Frequency range transmitter	768 - 872 MHz	f1 852 - 872 MHz	f1 818 - 872 MHz	768 - 872 MHz	768 - 785 MHz	852 - 872 MHz	840 - 872 MHz			
Frequency range receiver		f2 852 - 872 MHz	f2 768 - 818 MHz		852 - 872 MHz	w 2,4000 - 2,4835 GHz		818 - 872 MHz	768 - 872 MHz	850-872 MHz
Operating time	10 h	10 h	6 h	10 h	8 h	10 h	9 h	6 h	10 h	
Charging time	5 h	5 h	2 h	2 h	2 h	3 h	3 h	2 h	2 h	
Size	115 x 62 x 14 mm	115 x 62 x 14 mm	58 x 27 x 14 mm	84 x 39 x 14 mm	84 x 39 x 14 mm	97 x 41 x 17 mm	97 x 41 x 17 mm	60 x 51 x 20 mm	74 x 39 x 13 mm	8,5 x 9,2 x 12,6 mm
Weight	90 g	90 g	16 g	45 g	45 g	54 g	55 g	15 g	35 g	1,0 g
Range	30 m	30 m	30 m	30 m	30 m	30 m	30 m	30 m	30 m	25 m
Battery	Li-Polymer	Li-Polymer	Li-Polymer	Li-Polymer	Li-Polymer	Li-Polymer	Li-Polymer	Li-Polymer	Li-Polymer	Li-Polymer

Assortments may vary by land

### Recommended for:

Mobile phone	◆	◆					◆			
Landline telephone		◆								
Dialogue between 2 persons			◆	◆			◆			
Meetings and group work	◆	◆								
In the car			◆	◆	◆		◆			
Watching TV	◆	◆		◆	◆	◆	◆			
Integrated pupil (1 teacher)			◆	◆						
Integrated pupil (2 teachers)			◆	◆	◆					
Special school	◆	◆		◆	◆	◆	◆			
Lecture			◆	◆	◆					
Pupils with concentration problems				◆	◆					

The receiver is attached directly at the ear; can be used with or without all-in-the-ear hearing instruments.

Is used with hearing instruments, cochlear-implants, bone-anchored implants, head- or ear-phones; with induction loop.

Mini-receiver that can be used with hearing instruments, cochlear-implants and bone-anchored implants.